

PBMC Cryopreservation Media

CellStore

2026.1



CS-PM-D1 (10% DMSO, 100mL/bottle)

For peripheral blood– and cord blood–derived mononuclear cells and buffy coats.

Product Advantages

- Ready-to-use
- cGMP-manufactured with pharmacopeia ingredients
- Chemically Defined
- Animal component-free
- High-density, high-viability




Recommended Cell Density and Post-Thaw Performance

Density	Post-Thaw Viability		Viable Cell Recovery Rate
	3 Days	6 Months	
1 × 10 ⁶ – 1 × 10 ⁸ cells/mL	> 90%	> 90%	> 85%




Post-Thaw Viability and Recovery After Long-Term Cryostorage of PBMCs


Post-Thaw Viability



Viable Cell Recovery Rate



Post-Thaw Viability After 6-Month Liquid Nitrogen Storage




Preservation Density: 5×10^7 cells/mL; PBMCs from seven donors across multiple lots.

Note: Brand B refers to CryoStor® CS10 (BioLife Solutions).


- ❑ CS-PM-D1 delivers post-thaw cell viability comparable to fresh cells and higher than Brand B, with improved live cell recovery.
- ❑ These performance advantages are maintained over six months of liquid nitrogen storage.



Flow Cytometry Variability Often Starts before Analysis



Post-thaw Immune Cell Subset Distribution




Preservation Density: 5×10^7 cells/mL; immunophenotype analyzed by flow cytometry.

Note: Brand B refers to CryoStor® CS10 (BioLife Solutions).

- After thawing, both PBMC population structure and immune subset proportions (CD3⁺, CD19⁺, CD56⁺) in the CS-PM-D1 group were comparable to Fresh samples and more consistent than Brand B.




Whole Blood Frozen Directly with CS-PM-D1 (1:1) / No PBMC Isolation before Freezing / Stable Post-thaw Immune Profiles



- CS-PM-D1 enables direct whole blood cryopreservation without PBMC isolation, while maintaining stable post-thaw immune cell profiles.



Immune Cell Proportions Remain Stable



Why This Matters for Flow Labs

- ✓ Eliminates PBMC isolation before freezing
- ✓ Reduces handling-induced variability
- ✓ Enables consistent flow readouts across sample volumes

Particularly suitable for decentralized sample collection and centralized flow cytometry analysis.

- After thawing, the relative proportions of T cells (CD3⁺), B cells (CD19⁺), and NK cells (CD3⁻CD56⁺) remain consistent across storage volumes, supporting reliable downstream flow cytometry analysis.

CS-PM-D1 has been validated in flow cytometry–based workflows across research labs, biobanking, and clinical research settings.

Thanks